

[Remarkably
Everyday]

McGRATH[®]
MAC

Enhanced
Direct Laryngoscope



Distributed by



Maximize your chance of first-attempt success¹

Before you intubate, you prepare for every eventuality. So, whatever happens when your patient is on the table, when those critical moments are passing, you are better prepared.

Switch to the McGRATH[®] MAC enhanced direct laryngoscope (EDL) as your primary laryngoscope for every intubation. With one line of sight, you get direct and indirect visualization using the traditional direct laryngoscopy technique. With its lightweight, cord-free design and Macintosh blade, you get a laryngoscope you can use every day.

Make the move to the McGRATH[®] MAC enhanced direct laryngoscope as your primary laryngoscope, so when your patient is on the table, you are prepared.



McGRATH[®] MAC Enhanced Direct Laryngoscope



INTUITIVE

Uses standard direct laryngoscopy technique, while an anterior camera provides an enhanced view along a single line of sight.

PATIENT-FRIENDLY

Offers improved access to the patient's airway with a slim blade profile and without the excessive curve of most video laryngoscopes.²

CONFIDENCE

Includes a 250-minute battery with a minute-by-minute on-screen charge countdown that confirms device battery life.

DURABLE

Drop tested to two meters, the McGRATH[®] MAC EDL can endure rough handling in a fast-paced medical environment and still perform.

EASY TO CLEAN

Provides sterile integrity with disposable blades and is fully immersible. Both screen and handle can be treated with high-level disinfection.





ORDERING INFORMATION

McGRATH® MAC ENHANCED DIRECT LARYNGOSCOPE			
SKU NO.	DESCRIPTION	SIZE	QUANTITY
300-000-000	Handle	NA	1
350-017-000	McGRATH® MAC blade size 2	2	50/box
350-005-000	McGRATH® MAC blade size 3	3	50/box
350-013-000	McGRATH® MAC blade size 4	4	50/box
X3-003-000	McGRATH® MAC X blade size 3	3	10/box
340-000-000	250-minute battery pack	NA	1

TECHNICAL SPECIFICATIONS

McGRATH® MAC Handle

Size: 180 mm x 68 mm by 110 mm

Weight: 200g

Power: Proprietary 3.6V lithium battery pack (c. 250 minutes)

Light source: High-intensity LED

Display: 2.5" LCD color display

Camera: CMOS

Materials: Durable medical-grade thermoplastics with reinforced structural alloy core. The device and packaging are latex free.

McGRATH® MAC Disposable Blades

Material: Fog-free medical-grade optical polymer

Packaging: Packaged sterile for single use

For demonstrations and quotes, please contact your local representative. Details can be found at www.Covidien.com/YourNewPlanA19

McGRATH® MAC


 Designed & manufactured by Aircraft Medical

1. Su YC, Chen CC, Lee YK, Lee JY, Lin KJ. Comparison of video laryngoscopes with direct laryngoscopy for tracheal intubation: a meta-analysis of randomised trials. *Eur J Anaesthesiol.* 2011;28(11):788-795.

2. Levitan R, Heitz J, Sweeney M, Cooper R. The complexities of tracheal intubation with direct laryngoscopy and alternative intubation devices. *Ann Emerg Med.* 2011;57 (3):240-247.

Any comparative claims in this document are in respect to direct laryngoscopy.

Covidien is the sole acute care distributor for the McGRATH® MAC enhanced direct laryngoscope in the U.S., UK, Japan, Latin America, Australia, and New Zealand.

The McGRATH® MAC enhanced direct laryngoscope complies with EN 60601-1 and EN 60601-1-2 safety standards. The CE mark indicates that it meets the requirements of European Council Directive 93/42/ EEC concerning medical devices. The device is regulated in the USA under FDA Regulation Number 868.5540 and device listed under the name "McGRATH MAC".

"McGRATH" and "Aircraft" are registered trademarks of Aircraft Medical Limited. "CameraStick" is a trademark of Aircraft Medical Limited.

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. ©2013, 14 Covidien. 13-AW-0019(2)

Distributed by



6135 GUNBARREL AVENUE
BOULDER, CO 80301
800-635-5267

COVIDIEN.COM/RMS