

PRO 28Z™

Chair Cot



The **PRO 28Z Chair Cot** — provides exceptional maneuverability and stability. With its specialized ZFrame™ Advanced Positioning, the PRO 28Z transforms from a cot to a chair so that you can easily navigate tight hallways, narrow staircases, and elevators — all without transferring the patient. The PRO 28Z is compatible with the FERNO iNLINE® Cot Fastening System for SAE-compliant safety in the ambulance.

68% more STABILITY*



TIPRESIST™
ADVANCED STABILITY TECHNOLOGY

more MANEUVERABILITY



ZFRAME™
ADVANCED POSITIONING

- enables the PRO 28Z to transform from a cot to a chair for ultimate maneuverability in confined spaces
- **SAE J3027 Compliant** with iNLINE Cot Fastening System for in-ambulance safety
- **One-Hand release mechanism** quickly and easily switches the Pro 28Z between three locked positions: Chair, intermediate, and cot



Swing-up undercarriage

allows easy roll-in and roll-out unloading



Easily transitions from cot to chair position

with ZFrame Advanced Positioning



Cross-Compatible

When paired with the iNLINE Cot Compatibility Kit, the Pro 28Z is compatible with the SAE-compliant iNLINE Cot Fastening System

**Certified
SAE J3027
Compliant**

26g



**Cross-Compatible
System**

* Than the previous model, the 28Z PROFlexx, in the seated position.

PRO 28Z[™] Chair Cot



FEATURES AND BENEFITS

ENHANCED PATIENT CARE, COMFORT AND SAFETY

- **ZFrame Advanced Positioning** enables the Pro 28Z to transform from a cot to a chair, for ultimate maneuverability in confined spaces.
- **Most compact chair cot on the market** with a footprint of only 47" x 25" in the chair position.
- **One-hand release mechanism** quickly and easily switches the Pro 28Z between three locked positions: Chair, intermediate and cot.
- **TipRESIST Advanced Stability technology** provides 68% more tip resistance in the seated position*
- **Swing-up undercarriage** allows easy roll-in loading and roll-out unloading (One-man loading and unloading is possible when the cot is empty, two operators are always required during patient transports)

- **SX Surface Xtenders[™]** (Optional accessory) increase the patient surface area by 95%, giving a greater feeling of comfort and stability.
- **SAE-compliant Ferno X-Restraints** are color coded for ease of use

BUILT FOR LONGEVITY AND SIMPLE OPERATION

- Tested in a controlled environment to 10+ years of service
- **All-Metal Frame Construction** provides strength, durability, and performance with heavy patients.
- **2-Year Warranty** (Ask your Ferno Representative about our extended service agreements)

CROSS-COMPATIBLE DESIGN

- **iINLINE Compatibility Kit** allows you to secure the Pro 28Z to the iINLINE Cot Fastening System, making it SAE compliant. The Pro 28Z is also compatible with the 173 Cot Fastening System

SPECIFICATIONS

Item Number	Base model: 0012826
Length (folded)	81 in / 2063 mm
Length (raised)	76 in / 1925 mm
Width (overall)	25 in / 643 mm
Width (patient surface)	18 in / 448 mm
Height (folded)	8 in / 203 mm
Height	54.5 in / 1384 mm
Depth	47 in / 1202 mm
Weight**	96 lb / 43.5 kg
Weight Capacity	700 lb / 318 kg
Fastening Systems	Ferno iINLINE Cot Fastening System** Ferno 173 Cot Fastening System

STANDARD FEATURES

- Swing-Down Sidearms
- Ferno X-Restraints
- Two-Piece Locking Mattress
- Backrest Equipment Hook

COMPLIANCE***

The PRO 28Z has been designed to meet international test standards to increase safety for you and your patient.

- SAE J3027
- BS EN1789
- CAAS-GVS
- KKK-A-1822F
- NFPA 1917
- AMD 004
- IP56
- CE
- EN-1865



**Cot is weighed without mattress and accessories

***The Pro 28Z must be used with iINLINE[®] Cot Fastening System to be SAE compliant. When used with the Ferno 173 Cot Fastening System (ratchet-style fastener) it is not SAE compliant but IS AMD-compliant.